

BiliPOD 360 LED Intensive Care

Jaundice Management is critical for Newborns and premature babies. BILIPOD 360 offers a new MULTIDIRECTIONAL INTENSIVE PHOTOTHERAPY which brings new standards and benefits for newborns.



Ergonomic Design

Multi-direction Light Distribution

Complete exposure to light with Intensive Phototherapy care

Reduction in Treatment Time

Stable environment maintained along with temperature regulation with **integrated fan cooling system**



Electrical Elevation of the upper couple with the help of Up and Down keys present at the top of Graphic LCD Display

Slide out tray for easy access of the baby

High Irradiance Levels allow faster Bilirubin breakdown

Irradiance >60 uW/cm²/nm

Wavelength 450 – 480 nanometers

Lamp Unit Size: 70cm x 76cm x 165cm

Height Adjustability - up to 25cm

Weight - 50Kgs


Life of Lamps - 30000 hours

Graphic LCD Display allows display of Total Time, Baby Time and monitors Baby and Air Temperature



Manufactured by:

* AVI Healthcare Pvt Ltd reserves right to make changes without notice in design, specifications and models


AVI Healthcare Pvt Ltd
Manufacturing Neonatal Equipments "Glocally"

ISO 13485:2003
CERTIFIED COMPANY

ADDRESS: 25, Nanddeep Industrial Estate, Kondivita Lane,
Andheri East, Mumbai - 400059, Maharashtra, India
Telephone: +91 22 28320452/ +91 22 28326240
Email : sales@avihealthcare.com