



ISO 13485:2003 Certified

APNEA MONITOR

Neonatal Respiration Monitoring



AM 100 Apnea Monitor is capable of continuously monitoring the respiration of infants upto the age of 18 months. It gives audible warning and visual alarms if the infant stops breathing for more than preselected time.

Regular Breathing is indicated by Buzzer beep.

The Respiration Sensor consists of an encapsulated foam pad connected to a plastic tube. Expansion of the baby's abdomen (with each breath) will be detected by the unit.



TECHNICAL SPECIFICATIONS

DIMENSION	150 X 80 X 30 mm
WEIGHT	500gm without Adapter
POWER	9V (Rechargeable)
POWER ADAPTER	12V / 1A Battery
CONSUMPTION	10 Hours



Manufacturing Neonatal Equipments "Glocally"