



ISO 13485:2003 Certified

# BiliProbe Irradiance Meter

AVI Healthcare's BiliProbe Irradiance meter is designed to effectively measure accurate exposure level of standard neonatal phototherapy systems directly in the patient plane to produce optimal phototherapy results.

It is made for routine use in neonatal wards & nursing homes.



**Handy & Compact**

**Easy Operation**

**Bright LCD display**

**Accurate & Quick results**

## Technical Specifications

**Irradiance Range:** 0 - 200  $\mu\text{W}/\text{cm}^2/\text{nm}$

**Bandpass Filter Range:** 425 - 475 nm

**Measuring Time:** Approx. 4-5 seconds

**Minimal Graduation:** 1  $\mu\text{W}/\text{cm}^2/\text{nm}$

**Battery:** 9V Alkaline Battery

**Display:** LCD Screen Display

**Accuracy:**  $\pm 10\%$

**Press and Hold switch**

Manufactured by:

\* AVI Healthcare Pvt Ltd reserves right to make changes without notice in design, specifications and models

  
**AVI Healthcare Pvt Ltd**  
Manufacturing Neonatal Equipments "Glocally"

**ISO 13485:2003**  
CERTIFIED COMPANY

ADDRESS: 25, Nanddeep Industrial Estate, Kondivita Lane,  
Andheri East, Mumbai - 400059, Maharashtra, India

Telephone: +91 22 28320452/ +91 22 28326240

Email : sales@avihealthcare.com