



ISO 13485:2003 Certified

BiliProbe Irradiance Meter

AVI Healthcare's BiliProbe Irradiance meter is designed to effectively measure accurate exposure level of standard neonatal phototherapy systems directly in the patient plane to produce optimal phototherapy results.

It is made for routine use in neonatal wards & nursing homes.



Handy & Compact

Easy Operation

Bright LCD display

Accurate & Quick results

Technical Specifications

Irradiance Range: 0 - 200 $\mu\text{W}/\text{cm}^2/\text{nm}$

Bandpass Filter Range: 425 - 475 nm

Measuring Time: Approx. 4-5 seconds

Minimal Graduation: 1 $\mu\text{W}/\text{cm}^2/\text{nm}$

Battery: 9V Alkaline Battery


Display: LCD Screen Display

Accuracy: $\pm 10\%$

Press and Hold switch

Manufactured by:

* AVI Healthcare Pvt Ltd reserves right to make changes without notice in design, specifications and models


AVI Healthcare Pvt Ltd
Manufacturing Neonatal Equipments "Glocally"

ISO 13485:2003
CERTIFIED COMPANY

ADDRESS: 25, Nanddeep Industrial Estate, Kondivita Lane,
Andheri East, Mumbai - 400059, Maharashtra, India
Telephone: +91 22 28320452/ +91 22 28326240
Email : sales@avihealthcare.com