



ISO 13485:2003 Certified

# INFANTOMETER

AVI Healthcare's Infantometer is a infant height measurement device. It can be used everywhere and is particularly suited for use at doctors' clinics, hospitals as well as patients' home. It can be conveniently transported to any location

It has a flat firm surface board with fold up mechanism and low weight for babies upto 24 months.



## TECHNICAL SPECIFICATIONS

Measuring range: 45 cm - 90 cm

Graduation: 1 mm

Precision:  $\pm 0.2$ cm

Device weight: 1.2 kg

Unit of measurement: centimeters

Accuracy:  $\pm 0.2$ cm

Width of Board: 24 cm