



ISO 13485:2003 Certified

INFANTOMETER

AVI Healthcare's Infantometer is an height measurement device for infants up to 12 months. It can be used everywhere and is particularly suited for use at doctors' clinics, hospitals as well as patients' home. It can be conveniently transported to any location.

- ☀ It has thick, durable & heavy-duty ABS construction.
- ☀ It has flat, firm, sturdy, compact measurement surface with fold-up mechanism.
- ☀ It ensures precise & accurate measurement with easy to read double sided parallel positioned scale.
- ☀ No need for calibration or set-up as it has adjustable tight & firm connection links.



TECHNICAL SPECIFICATIONS

Measuring range: 45 cm - 100 cm

Graduation: 0.1 cm

Precision: ± 0.2 cm

Device weight: 1.2 kg

Unit of measurement: Centimeter

Accuracy: ± 0.2 cm

Width of Board: 30 cm (approx.)

Manufactured by:

* AVI Healthcare Pvt Ltd reserves right to make changes without notice in design, specifications and models


AVI Healthcare Pvt Ltd
Manufacturing Neonatal Equipments "Glocally"

ADDRESS: 25, Nanddeep Industrial Estate, Kondivita Lane,
Andheri East, Mumbai - 400059, Maharashtra, India
Telephone: +91 22 28320452/ +91 22 28326240
Email : sales@avihealthcare.com

ISO 13485:2003
CERTIFIED COMPANY