



BiliPOD Portable LED Double Surface

Uniform Light Distribution
with 11 High Intensity LED Lamps and Bilipad

Lamp Unit Size 45cm x 30cm x 5cm

Lamp Tilt up to 90°

Life of Lamps - 30000 hours


Zero Noise Cooling Fan

Effective Irradiance >30μW/cm²


Low Power Consumption

Easy Installation & Low Maintenance

Optimum Wavelength of 460nm



High Irradiance levels
allows faster Bilirubin Breakdown



BILIPAD LED Phototherapy is used for undersurface phototherapy

Light Weight
Lamp Unit weighs not more than 3 Kg.

Observation Lamp
Two Observation White LEDs provide enough light to observe the color of baby's skin



Modern Design

It is designed to occupy minimal space without compromising the comfort & treatment efficacy.



Portable and Compact

It can be carried / placed near mother's bed.

Takes minimum space and the lamp unit is very sleek


LCD based Timer

Total Lamp Usage - hhhhh
Patient Timer - hhh:mm



Manufactured by:

* AVI Healthcare Pvt Ltd reserves right to make changes without notice in design, specifications and models


AVI Healthcare Pvt Ltd
Manufacturing Neonatal Equipments "Glocally"

ISO 13485:2003
CERTIFIED COMPANY

ADDRESS: 25, Nanddeep Industrial Estate, Kondivita Lane,
Andheri East, Mumbai - 400059, Maharashtra, India

Telephone: +91 22 28320452/ +91 22 28326240

Email : sales@avihealthcare.com