



ISO 13485:2003 Certified

STADIOMETER

AVI Healthcare's Stadiometer is a mobile height measurement device. It can be used everywhere and is particularly suited for use at doctors' clinics, hospitals as well as patients' home. It can be conveniently transported to any location.



ACCURATE RESULTS

PORTABLE

EASY TO TRANSPORT

EASY TO USE

WELL DESIGNED

LIGHT WEIGHT

RELIABLE

SIMPLE ASSEMBLY

TECHNICAL SPECIFICATIONS

Measuring range: 20 – 205 cm / 8 – 81"

Graduation: 1 mm / 1/8"

Dimensions(WxHxD): 328 x 2,145 x 574 mm

Device weight: 3.6 kg