

STADIOMETER

AVI Healthcare's Stadiometer is a portable height measurement device. It can be used everywhere and is particularly suited for use at doctors' clinics, hospitals as well as patients' home. It can be conveniently transported to any location.

- ☀ It has firm vertical stand with easy fold-up for simple set-up with firm & tight fixations.
- ☀ It has a large skid-proof stable foot-plate with spirit level indicator for better judgement with blunt edges/corners for safer usage.

Salient Features

Skid-proof Foot plate

Smooth measuring wedge

Spirit-level indicator

Durable & Heavy duty

Resistant to heat & humidity

Simple Assembly

Portable

Accurate Results

Easy to Transport

Easy to Use

Well Designed

Reliable

TECHNICAL SPECIFICATIONS

Measuring range: 20 – 205 cm / 8 – 81"

Board depth: 5.5 cm

Graduation: 0.1 cm

Precision: ± 0.2 cm

Unit of Measurement: Centimetre

Accuracy: ± 0.2 cm

Dimensions(WxHxD): 32.8 x 214.5 x 57.4 cm

Device weight: 3.6 kg

* AVI Healthcare Pvt Ltd reserves right to make changes without notice in design, specifications and models

Manufactured by:

ISO 13485:2016
CERTIFIED COMPANY



Manufacturing Neonatal Equipments "Glocally"

Address : 25, Nanddeep Industrial Estate, Kondivita Lane,
Andheri East, Mumbai - 400059, Maharashtra, India

Telephone : +91 22 28320452/ +91 22 28326240

Email : sales@avihealthcare.com

