

# AVIA TransVent

## Neonatal Transport Ventilator

AVIA TransVent - Neonatal Transport Ventilator is a portable, light weight, robust, and is able to function in demanding environments with little maintenance.

When the smallest lives are in your hands, you can depend on AVIA TransVent emergency and transport Pediatric / Neonatal Ventilator.

**TRANSPORT**

**COMPACT**

**PORTABLE**

**DEPENDABLE**



**EASY TO USE**

**MODES:  
CMV, CPAP**

**CALIBRATED  
CONTROLS**

**LIGHT WEIGHT**

### Technical Specifications

<b>OPERATION</b>	Time Cycled, Pressure Limited, Continuous Flow Ventilator	<b>PRESSURE SCALE</b>	-20 to 100 cmH <sub>2</sub> O
<b>MODES</b>	CPAP, CMV	<b>AIR FLOW</b>	0 to 15 lpm
<b>INSPIRATION TIME</b>	0 to 2 sec	<b>OXYGEN FLOW</b>	0 to 15 lpm
<b>EXPIRATION TIME</b>	0 to 6 sec	<b>TIDAL VOLUME</b>	330 ml to 0
<b>PEEP</b>	0 to 50 cmH <sub>2</sub> O	<b>BATTERY BACKUP</b>	Upto 4 hours
<b>PIP</b>	0 to 30 cmH <sub>2</sub> O	<b>DIMENSION</b>	270 mm x 195 mm x 150 mm
		<b>WEIGHT</b>	5 Kgs